



Sterifast

Sterilization & Disinfection

STERILIZATION PACKAGING

STERILIZATION PACKAGING

The end user of our product divided into two groups:

1. Hospital, Clinics, and healthcare field,...

*According to EN 868 -3, EN 868 -5 and ISO 11607

2. Medical Device Manufacturers

*According to EN 886 -6, EN 686 -5 and ISO 11607



OWN BRAND LABELING (OBL) :

Sterifast offers medical device manufacturers their packages with favorite design and brands.

Printing on medical packages includes, Indicators (Steam, E.O., VH_2O_2), Customers'

designs and Customers' Information (product' s name, logo, technical specifications, standards and etc.)



TYVEK STERILIZATION ROLL

- Color: Transparent
- Indicator: VH_2O_2
- Customized size



FLAT STERILIZATION ROLL AND POUCH

- Color: Blue , Green , Transparent
- Indicators: Steam & E.O. and Formaldehyde upon request
- Customized size



MEDICAL FILM

- Different thicknesses
- Color: Transparent
- Customized solution based on PA, PE, PP,...
- Customized width, Customized length



GUSSETED STERILIZATION ROLL

- Color: Blue, Green, Transparent
- Indicators: Steam, E.O. and Formaldehyde upon request



SELF-SEALING POUCH

- Color: Blue, Green
- Customized size
- Indicators: Steam & E.O. and Formaldehyde upon request



ALUMINIUM POUCH AND ROLL

- Different layers with combination of paper, PE & AL
- Customized width, Customized length



TYVEK[®] HEADER BAG

Tyvek header bags consist of see-through or white film and a Tyvek lap. This allows complete visibility of these products contained inside and makes pouches ideal for Custom-Pack suitable for E.O. and irradiation sterilization.
 - Indicators: E.O. and irradiation upon request.



SMS SHEET

SMS non woven:
 - 100% synthetic fibers
 - Minimum of 3 layers of Spun bond/ Melt blown / spun bond



CREPE PAPER

- High Quality Creped Paper For Sterilization Packaging
 - Standard crepe: 100% cellulose material
 - Medical crepe paper works as an efficient barrier to humidity and it can be used to pack
 - Materials that will undergo either steam or E.O sterilization



LABEL GUN & INDICATOR LABEL

- Conform to ISO11140 -1 (Type 1)
 - Indicator: steam or VH_2O_2
 - Printing operator code, cycle number, expiry date, etc. by label gun



BIOLOGICAL INDICATOR

Biological indicators are test systems that contain viable microorganisms with a defined resistance to a specific sterilization process. They help monitor whether the necessary conditions were met to kill a specified number of microorganisms for a given sterilization process.



AUTOMATIC SEALING MACHINE

- Automatic Sealing Machine
 - With an integrated printer
 - Quick, safe and traceable sealing

STERIFAST CUTTER MACHINE

- Ergonomic design
 - Useful length: 50, 70, 90 (cm)



TYPE 4 (14)
CHEMICAL INDICATOR



TYPE 6 (16)
CHEMICAL INDICATOR



EO INDICATOR (E1)



BOWIE & DICK TEST



Sterifast

Sterilization & Disinfection

ecoPlasma